

**MSC-21947-1-CU**

**Express mailing certificate No. EU657647915US**  
**Patent Application**

**MOLECULAR-SPECIFIC UROKINASE ANTIBODIES**

**ORIGIN OF THE INVENTION**

**[0001]**        The invention described herein was made at least in part in the performance of work under a NASA contract and is subject to Public Law 96-517 (35 U.S.C. §200 et seq.). The contractor has not elected to retain title to the invention.

**Field of the Invention**

[0002] The invention relates to compositions and methods of using antibodies to analyze molecular forms of urokinase. More particularly, a group of such antibodies, each of which recognizes and selectively binds to one of several molecular forms of urokinase, is used to detect, quantitate or purify individual forms of urokinase. Species of urokinase that can be analyzed using the compositions and methods of the invention include the inactive 54,000 dalton single chain form, the bioactive, 34,000 dalton double-chain form and the bioactive 33,000 dalton form.

**Description of the Related Art**

- [0003] Urokinase is one member of a class of proteins known as fibrinolytic proteins. Other fibrinolytic proteins include streptokinase and tissue plasminogen activator (tPA). These type proteins aid in natural blood clot dissolution. Blood clot dissolution is an important therapy in the treatment of several major debilitating and life threatening diseases including myocardial infarction associated with heart attacks, stroke, deep vein thrombosis and pulmonary embolus.
- [0004] Urokinase is more specifically a serine protease. Its natural substrate is plasminogen, which it converts to bioactive plasmin. When plasmin comes in contact with the major structural component of blood clots called fibrin, the fibrin is degraded into small fragments dissolving the clot.
- [0005] Consequently, there is great interest in urokinase as a therapeutic agent. U.S. Patent No. 4,873,083 discloses a fibrinolytic composition employing an enzyme such as urokinase tissue plasminogen activator and a surface active copolymer for the treatment of thrombosis. U.S. Patent Nos. 4,600,580, 4,673,573, and 4,741,903 disclose chemically-modified forms of urokinase that employ urokinase for fibrinolytic activity in therapeutic applications. U.S. Patent No. 4,791,068 discloses the use of a polyclonal antibody system for detecting an inhibitor of urokinase-type plasminogen activator and its use in diagnostic assays.
- [0006] For these reasons, there is also great interest in developing more efficient methods of urokinase production. U.S. Patent No. 4,889,808 discloses a method of enhancing the production of a single chain urokinase by the addition of heparin and a cell growth factor to culture medium of cells producing the enzyme. The production levels were determined using polyclonal goat antisera which was not specific for active forms of urokinase, but did allow crude detection of enhanced levels of production.

- [0007] Another approach is to isolate novel cell lines that produce large quantities of this protein. Urokinase is known to be secreted by kidney cells, some neoplastic tissues, as well as embryonic cell line (Bernick and Kwann, *J. Clin. Invest.* 48:1740-1753 (1969); Barlow et al., *Thrombosis Res.*, 32:29-34 (1983)). Bernick and Kwann (*J. Lab. Clin. Med.* 70:650 (1967)) have reported that only 5-10% of the cells from human embryonic kidney (HEK) produce plasminogen activators. Enhanced production can be achieved if cells producing urokinase are isolated from those cells that do not produce urokinase. When a heterogeneous sample of renal cells is separated under microgravity conditions by continuous flow electrophoresis, several distinct subgroups can be isolated. To determine which subgroup most efficiently produces urokinase, sensitive and specific assay methods that can detect the various molecular forms and total urokinase are required.
- [0008] Since urokinase is a cancer-specific protein, it is a marker of neoplastic cell development, particularly in cases of invasion and metastasis where this enzyme can play an important role. The role of urokinase in cancer has been discussed as it relates to breast cancer (Schmitt et al., *Blood Coagulation and Fibrin.* 1:695 (1990), urinary cancer (Hasui et al., *Cancer Res.* 49:1067 (1989), and prostate cancer (Gaylis et al. *J. Urology* 142:193 (1989)). Thus, sensitive and specific assays are also needed which can be used to diagnose these conditions.
- [0009] Urokinase is synthesized as a single polypeptide of 431 amino acid residues. An amino terminal signal peptide of 20 amino acids (-20 to -1) is cleaved from the polypeptide as it is processed during export from the cell, leaving a 411 amino acid protein (1-411) that is called the urokinase zymogen or single chain urokinase plasminogen activator (scuPA).
- [0010] Urokinase zymogen has very little or no enzymatic activity. Activation of the urokinase zymogen involves cleavage of the peptide bond between lysine 158

and isoleucine 159 and the loss of lysine 158 to form the high molecular weight (HMW) 54,000 dalton active enzyme (Verstraete and Collen, *Blood* 67:1529-1541 (1986). Activation can also include (in addition to 158-159 cleavage) cleavage of the bond between lysine 135 and lysine 136 followed by loss of the amino terminal 135 amino acids (the "A chain") leaving an active low molecular weight (33,000 dalton) enzyme. The non-enzymatic products of these cleavage events are rapidly degraded. Activation and autodegradation of urokinase in a cell culture can result in all three types of urokinase being present in culture medium at one time.

[0011] To identify cell populations that express increased levels of urokinase a quantitative measurements of all forms of urokinase is necessary. To be of practical use, such assays must be sensitive enough that urokinase production from small numbers of cells can be analyzed. Such assays must also be relatively easy to perform so that large numbers of samples containing as many as  $5-10 \times 10^4$  cells/ml or more can be tested.

[0012] Current urokinase assay methods are cumbersome. One of the common procedures is the fibrin plate lysis assay. See, e.g., U.S. Patent 4,741,903; and, Lewis et al., "Problems in the Bioassay of Products from Cultured HEK Cells: Plasminogen Activator" in *Eukaryotic Cell Cultures: Basics and Applications* 172 Plenum Press, New York 1984. This assay measures fibrinolytic activity. Test samples are incubated with fibrin clots and the increase in the diameter of the clot lysis zone is measured as a function of time. The rate of clot disappearance is a measure of fibrinolytic activity. The detection limit for such assays is generally 10-15 I.U./ml. Other assays include the timed clot-lysis assay (Swank et al. *Computer Meth. Programs Biomed.*, 33:95-105 (1990)),  $^{125}\text{I}$ -labeled fibrinogen assay and colorimetric assays which use synthetic chemical substrates, and the S-2444 assay, KabiVitrum, Sweden (Claeson et al., *Haemostasis* 7:76-78 (1978). Assays such as these are also activity assays since they are based on measurement

of fibrinolytic activity. They are used to detect and quantitate active urokinase but are incapable of detecting the inactive zymogen.

[0013] Furthermore, existing assays are inadequate to detect small quantities of uPA in very small (e.g. 10  $\mu$ l) samples of culture medium or bodily fluids. Thus, such assays may not identify potentially important cell lines that produce large quantities of the urokinase zymogen (scuPA). The scuPA levels can be approximated by preincubation of samples with plasmin (which activates urokinase zymogen) followed by an activity assay. However, additional steps in an already cumbersome assay increase the variation in assay results. Moreover, it is difficult to determine the amount of plasmin and hydrolysis times needed to fully activate zymogen in all cases.

[0014] Inconsistent test results are frequently associated with fibrin clot activity assays. For example, storage of samples affects apparent activity, dose response curves are nonlinear, and replicate cell cultures produce variable activities. More importantly, protein diffusion is a factor affecting the apparent activity of urokinase. The 33,000 dalton form diffuses faster than the 54,000 dalton form increasing the apparent activity relative to the 54,000 dalton form. Although the chromogenic assay is not limited by protein diffusion, this assay is not as sensitive as the fibrin clot assay. Sample composition can greatly affect activity measurements, for example the presence of inhibitors of urokinase may greatly reduce the amount of activity measured for a given amount of enzyme. Similarly, other non-urokinase fibrinolytic enzymes in the sample may interfere by increasing the apparent activity. Moreover, these assays are not sensitive enough to detect urokinase produced from different subpopulations of small numbers of cells.

[0015] Activity assays also are unable to distinguish between various active forms of urokinase since they only measure total fibrinolytic activity. This makes comparisons between samples difficult. Use of a combination of assays, such as

the clot lysis assay in combination with the chromogenic assay, can give approximate estimates of the ratios of the 54,000 and the 33,000 dalton forms of urokinase. However, this is a cumbersome and time-consuming procedure.

[0016] An alternative approach to analysis of the urokinase from cell culture is to use quantitative immunological techniques such as enzyme linked immunosorbent assays (ELISA) or radioimmunoassays (RIA). Various murine and rabbit IgG monoclonal antibodies are available against the A-chain (residues 1-135) or B-chain (residues 159-411) of urokinase (i.e., American Diagnostica, New York). These antibodies are made from animals immunized with the major molecular forms of urokinase. U.S. Patent No. 4,474,756 discloses a method of producing antibodies by immunization of animals with a whole protein, such as human urokinase, after chemically crosslinking the protein to a polysaccharide. However, since urokinase has more than 70% homology with trypsin and significant homology to other serine proteases and tissue plasminogen activator (tPA), specificity is a problem. Therefore, non-specific antibodies such as these also bind to and detect other serine proteases and fibrinolytic proteins, and they do not give reliable quantitative estimates of urokinase in samples.

[0017] Thus, compositions and methods are needed that can sensitively and specifically quantitate urokinase in all its forms in a sample. Assays are needed which overcome the cumbersome aspects of existing methods and which produce less variable results and which are not subject to interferences from inhibitors which could contaminate sample preparations. Assays are needed which are specific for urokinase and which distinguish between closely related proteins such as tissue plasminogen activator or streptokinase or homologous serine proteases that increase the "apparent" concentration of urokinase observed in samples making existing assays inaccurate. Lastly, assays are needed which are capable of measuring the concentration of each form of urokinase in a sample.

**SUMMARY OF THE INVENTION**

[0018] In accordance with the present invention, polyclonal and monoclonal antibody compositions are provided that can sensitively and specifically quantitate urokinase in all its forms. The compositions provided can be used in methods that are convenient to perform, give reliable results, and are free from interferences from enzyme inhibitors that could affect fibrinolytic activity measurements. The compositions are specific for urokinase and will not detect closely related proteins such as tissue plasminogen activator or homologous serine proteases. Compositions provided with this invention can be used to measure the concentration of each form of urokinase even in a small sample.

[0019] Antibodies (polyclonal and monoclonal) were generated against urokinase by employing synthetic peptides as immunogens. The antibodies that are generated by injecting an immunogen into an animal are said to be "directed against" that immunogen. Antibodies directed against immunogens may bind to other molecules as well. Particularly, antibodies directed against an immunogen may bind to other biological macromolecules that contain the immunogen. In particular, where as here the immunogen is a peptide which represents a particular portion of the urokinase macromolecule, antibodies directed against the peptide immunogen can bind to the biological macromolecule urokinase in any of its forms which contains the immunogen peptide sequence. If the immunogen is found on one macromolecule or a subset of macromolecules, and the antibody directed against that immunogen binds preferentially to that macromolecule or subset of macromolecules, the antibody is said to have "specific" binding for that macromolecule or subset of macromolecules. "Specific" binding is merely a reflection of the fact that the binding affinity of the antibody to macromolecules



containing the immunogen is higher than the binding affinity of the same antibody to other molecules.

**[0020]** In certain aspects, the invention relates to an antibody directed against a peptide whose sequence is chiefly derived from the urokinase amino acid sequence (Seq. ID No. 16). "Chiefly derived" for purposes of this invention means that the sequence of the peptide is virtually identical to at least a portion of a linear sequence of urokinase except for minor alteration such as additional terminal residues or substitutions of cysteine residues to prevent internal disulfide linkages (i.e., see Seq. ID No. 5 where cysteine-131 is replaced by a glycine residue in the peptide to prevent disulfide formation with cysteine-126).

**[0021]** The antibody is first characterized by its capability to specifically bind to urokinase. As further amplified in the claims, the characteristic of specific binding to urokinase is such that the antibody possesses a substantially lower binding affinity to proteins other than urokinase than that it has for urokinase. Importantly, this reduced binding affinity for non-urokinase proteins includes those proteins as similar in amino acid sequence to that of urokinase as is trypsin, a protein with a high degree of sequence homology with urokinase.

**[0022]** One group of antibodies that meet this general description will be antibodies which bind to the amino terminal end of urokinase. The amino terminal end of the urokinase molecule for purposes of the invention comprise urokinase amino acid residues 1-135 of Seq. ID No. 16.

**[0023]** Among those antibodies which are specific for the amino terminal end of the urokinase molecule will be those antibodies directed against one or more of six specific peptides. These peptides may be any one of a number of such peptides based upon sequences of the amino terminal sequence of urokinase such as: SNELHQVPSNCD (Seq. ID No. 1), RGKASTDTMGRPCLP (Seq. ID No. 2), CRNPDNRRRP (Seq. ID No. 3), RNPDNRRRPWC (Seq. ID No. 4),

CMVHDGADGK (Seq. ID No. 5) or MVHDCADGK (Seq. ID No. 6). In certain instances, an original residue found in urokinase will be altered for ease of synthesis which alteration does not materially affect the immunological response to the synthetic peptide (e.g., instead of asparagine-10 as in urokinase, Seq. ID No. 1 possesses an aspartate; instead of cysteine-131, Seq. ID No. 5 possesses a glycine).

[0024] Another group of antibodies that meet this general description will be antibodies which bind to the carboxyl terminal end of urokinase. The carboxyl terminal end of the urokinase molecule for purposes of the invention comprise urokinase amino acid residues 159-411 of Seq. ID No. 16.

[0025] Among those antibodies which are specific for the carboxy terminal end of the urokinase molecule will be those antibodies directed against one or more of six specific peptides. These peptides may be any one of a number of such peptides based upon sequences of the amino terminal sequence of urokinase such as: YRRHRGGSVTYVC (Seq. ID No. 7), CFIDYPPKKEDY (Seq. ID No. 8), SRLNSNTQGEMK (Seq. ID No. 9), SMYNDPQFGTSC (Seq. ID No. 10), LISHRECQQPHYYGSEVTTKMLC (Seq. ID No. 11), or SHTKEENGLAL (Seq. ID No. 12).

[0026] Other antibodies which meet the general descriptions provided by the antibodies of the invention include antibodies which can detect cleavage sites in urokinase. These antibodies include an antibody which has a binding site that includes the peptide bond between lysine 158 and isoleucine 159. More specifically, such an antibody will be one which is directed against one or more of two peptides selected from the group of peptides comprising PRFKIIG (Seq. ID No. 13) or RPRFKIIGGE (Seq. ID No. 14). These antibodies also include an antibody which has a binding site that includes the peptide bond between lysine 135 and lysine 136. More specifically, such an antibody will be one which has a

lower binding affinity for urokinase when no chemical bond exists between lysine 158 and isoleucine 159. Similarly, such an antibody will be one which has a lower binding affinity for urokinase when no chemical bond exists between lysine 135 and lysine 136. The site defined by lysine 135 and lysine 136 is the natural cleavage site of the urokinase proenzyme and is structurally exposed to cleavage, and thus is a likely natural epitope region. This site is bridged by Seq. ID No. 17.

[0027] Particularly, in a preferred embodiment, an antibody directed against the peptide sequence RPRFKIIGGE (Seq. ID No. 14), with a binding affinity constant for the urokinase zymogen of at least  $1 \times 10^8 \text{ M}^{-1}$ , is disclosed which additionally exhibits the binding affinity constant for urokinase lacking a peptide bond between lysine 158 and isoleucine 159 decreased by at least ten-fold.

[0028] It is also possible to make multimers of peptides by covalently linking one or more of peptides in Seq. ID Nos. 1-15. Thus, duplex peptides may easily be covalently linked by disulfide linkages, as for example by making (102-111)—(126-135), (59-73)—(103-113), (177-189)—(205-215). It is also possible to cause a peptide to form an internal disulfide, as for example in 319-341 (Seq. ID No. 11).

[0029] It will be recognized by one of skill in the art of making antibodies against specific peptides, that such antibodies may be derived from polyclonal sera. Alternatively, and in certain aspects, preferably, the antibody of any of the embodiments mentioned above may be a monoclonal antibody. More specifically, the monoclonal antibodies of the invention have been given the designations: A3-E5-A7, A3-E5-B8, A3-E5-G10, A3-E5-H12, A3-E5-B3, C2-A4-H4, C2-A4-F10, C2-A4-F11, D4-A5-E7, D4-C1-A5, D4-C1-G4, D4-C6-A4, D4-C6-A8, D4-C6-E12, D4-C6-F7 or D4-C6-G12. Hybridoma cell lines capable of producing such monoclonal antibodies are also provided. While these specific monoclonal antibodies and the hybridomas that produce them represent the results of the particular experiments defined in greater detail below, they are in merely

exemplary in nature. Using the sequences selected and disclosed in the invention, the normally skilled practitioner may readily obtain antibodies with the same characteristics as reported here and which will function identically to those reported here.

**[0030]** The peptides will be one of the group of peptides disclosed in Sequence ID Nos. 1-15. It will be recognized, however, by those of skill in the art that the catalysts described above and those claimed in general may contain functionally equivalent amino acid substitutions. The importance of the hydropathic index of amino acids in conferring biological function on a peptide is generally known by those of skill in the art. It has been found by many researchers that certain amino acids may be substituted for other amino acids having a similar hydropathic index or score and still retain similar if not identical biological activity. As displayed below, amino acids are assigned a hydropathic index on the basis of their hydrophobicity and charge characteristics. It is believed that the relative hydropathic character of the amino acid determines the secondary structure of the resultant peptide, which in turn defines the interaction of the peptide with the substrate molecule. It is proposed that biological functional equivalence may typically be maintained where amino acids are exchanged having no more than a  $\pm 1$  to 2 difference in the index value, and more preferably within a  $\pm 1$  difference.

<b>AMINO ACID</b>	<b>HYDROPATHIC INDEX</b>
Isoleucine	4.5
Valine	4.2
Leucine	3.8
Phenylalanine	2.8
Cysteine/Cystine	2.5
Methionine	1.9
Alanine	1.8
Glycine	-0.4
Threonine	-0.7
Tryptophan	-0.9
Serine	-0.8
Tyrosine	-1.3
Proline	-1.6
Histidine	-3.2
Glutamic Acid	-3.5
Glutamine	-3.5
Aspartic Acid	-3.5
Asparagine	-3.5
Lysine	-3.9
Arginine	-4.5

[0031] Thus, for example, isoleucine, which has a hydrophatic index of +4.5, can be substituted for valine (+ 4.2) or leucine (+ 3.8), and might still obtain a peptide

having similar biological activity. Alternatively, at the other end of the scale, lysine (-3.9) can be substituted for arginine (-4.5), and so on.

**[0032]** Accordingly, these amino acid substitutions are generally based on the relative similarity of R-group substituents, for example, in terms of size, electrophilic character, charge, and the like. In general, although these are not the only such substitutions, the preferred substitutions which take various of the foregoing characteristics into consideration include the following:

<b>ORIGINAL RESIDUE</b>	<b>EXEMPLARY SUBSTITUTIONS</b>
alanine	glycine; serine
arginine	lysine
asparagine	glutamine; histidine
aspartic acid	glutamic acid
cysteine	serine
glutamine	asparagine
glycine	alanine
histidine	asparagine; glutamine
isoleucine	leucine; valine
leucine	isoleucine; valine
lysine	arginine; glutamine; glutamic acid
methionine	leucine; tyrosine
serine	threonine
threonine	serine
tryptophan	tyrosine
tyrosine	tryptophan; phenylalanine
valine	isoleucine; leucine

**[0033]** In addition to the compositions provided by the invention, methods are disclosed which utilize the anti-urokinase antibody compositions in a variety of ways. Thus, in one embodiment, a method of detecting urokinase zymogen in a sample is provided. This method comprises standard plate assays including solid phase plate radioimmunoassay or enzyme linked immunoassays (ELISA).

Another method made possible by the anti-urokinase antibody compositions of the invention include detecting urokinase associated with the surface of a cancer cell. This method comprises using the antibodies of the invention to detect urokinase levels in a cell population, which levels are correlated to DNA content for evaluation of potential of tumor cells.

**[0034]** For whatever purpose the quantitation is accomplished, the technique will be similar in its use of the peptides of the invention. Thus, Seq. ID Nos. 1-15 can be used to generate immunological compositions capable of recognizing the urokinase zymogen. The same is true of a peptide comprising residues 135 and 136 of Seq. ID No. 16. When the same battery of immunological compositions is used against the high molecular form of urokinase which lacks the 158-159 peptide bond, each such composition will bind the high molecular weight form of urokinase except immunological composition derived from Seq. ID Nos. 13 and 14, thusly providing a test to distinguish the inactive zymogen from the active high molecular weight form of urokinase. The same battery of immunological compositions is useful in distinguishing the low molecular weight form will not efficiently bind immunological compositions derived from Seq. ID Nos. 1-6, 13-14 or from a peptide comprising residues 135 and 136 of Seq. ID No. 16.

**[0035]** Diagnostic kits for detecting the presence and quantity of urokinase in a sample are also provided. Such kits comprise effective amounts of the antibodies which individually bind to the amino terminal end of urokinase, the carboxyl terminal end of urokinase, the binding site that includes the peptide bond between lysine 158 and isoleucine 159, and the binding site that includes the peptide bond between lysine 135 and lysine 136. Additionally, each such antibody includes an indicating group capable of detecting binding of the antibody to a target molecule and standards of an effective concentration of each of the bioactive and inactive forms of urokinase detectable using the four antibodies included in the kit.



**BRIEF DESCRIPTION OF THE DRAWINGS**

[0036]       **Figure 1:**     Standard curve for Run 4 Day 12 supernatants radioimmunoassay.

[0037]       **Figure 2:**     Bar graph comparing amount of scuPA in the Run 4 Day 12 supernatants using two prior art methods (S-2444 and MCLA) and methods of invention (MSA-RIA).

**DESCRIPTION OF PREFERRED EMBODIMENTS**

[0038] With respect to the inactive form of urokinase (See Seq. ID No. 16), activation of the zymogen involves cleavage of the peptide bond 158-159 and the loss of Lys-158. An active form of the enzyme is produced which has the sequence 1-157 (A chain) linked to sequence 159-411 (B chain) by a disulfide bond between Cys-148 and Cys-249. Another active form of the enzyme has sequence 136-157 (A1 chain) linked to sequence 159-411 (B chain). The various regions selected for synthesis correspond to sequences 1-12 (Seq. ID No. 1), 59-73 (Seq. ID No. 2), 102-111 (Seq. ID No. 3), 103-113 (Seq. ID No. 4), 126-135 (G at 131) (Seq. ID No. 5), 127-135 (Seq. ID No. 6), 177-189 (Seq. ID No. 7), 205-215 (Seq. ID No. 8), 222-233 (Seq. ID No. 9), 319-330 (Seq. ID No. 10), 319-341 (Seq. ID No. 11), 401-411 (Seq. ID No. 12), 155-161 (Seq. ID No. 13) and 154-163 (Seq. ID No. 14). In addition, the following peptide pairs will be linked by disulfide bonds: (59-73), (102-111)-(126-135, G131) and (177-189)-(205-215). The covalent structures of the peptides that have been (or are being) synthesized are shown in Table 1. The single letter notations of the amino acids are: A, alanine; C, cysteine; D, aspartic acid; E, glutamic acid; F, phenylalanine; G, glycine; H, histidine; I, isoleucine; K, lysine; L, leucine, M, methionine; N, asparagine; P, proline; Q, glutamine; R, arginine; S, serine; T, threonine; V, valine; W, tryptophan; Y, tyrosine.

Peptide Immunogen Synthesis

- [0039] The peptides used for certain particular examples of the invention correspond to the Sequence ID Nos. 1-6. These peptides, as well as other peptides produced using the methods of the invention, are synthesized either by t-butyloxycarbonyl (t-Boc) on a phenylacetamidomethyl (PAM) resin or by 9-fluorenmethylcarbonyl (Fmoc) amino acids on a benzyloxymethyl alcohol resin as described elsewhere (McCormick and Atassi, *Biochem. J.*, 224:9950-10000 (1984); Mulac-Jericevic and Atassi, *J. Prot. Chem.*, 6:365-373 (1987); Atassi et al., *Proc. Natl. Acad. Sci. USA*, 88:3613 (1991), together with methods for purification and characterization of the peptides.
- [0040] Peptide immunogens can be synthesized chemically by the well known methods discussed above. They can also be isolated by purification after proteolysis of individual forms of urokinase. Particular care must be given to the choice of protease used in the latter method so that desired epitomes are not destroyed.
- [0041] To avoid generating antibodies that cross react with other known proteins peptide immunogen sequences are selected which are unique to urokinase. To determine which sequences are unique to urokinase, the urokinase sequence can be aligned with other proteins whose sequences are known. Many of these protein sequences can be found in data bases and convenient computer programs can be enlisted to perform a homology search. Such searches may be accomplished using the FAST-DB (Fast Pairwise Comparison of Sequences - Release 5.4) sequence homology searching program as provided by Intelligenetics, Inc., 700 East El Camino Real, Mountain View, CA 94040. The searches may be conducted in the publicly-available sequence databases -- PIR39 and Swiss-Prot 28.
- [0042] Fourteen peptides were prepared for use as immunogens. Peptides derived from the amino terminal end of the urokinase protein, including residues 1 through

135, were used to generate antibodies specific for the A chain of urokinase. Six peptides represent areas present only in the zymogen and high molecular weight forms of the enzyme. The low molecular weight form does not contain these areas due to the loss of the first 158 residues after enzymatic degradation. As a result, antibodies directed against these peptides only bound the high molecular weight forms of urokinase, providing compositions for simple discriminating assays. The six peptides that were found useful in generating the A chain specific antibodies include: SNELHQVPSNCD (Seq. ID No. 1), RGKASTDTMGRPCLP (Seq. ID No. 2), CRNPDNRRRP (Seq. ID No. 3), RNP DNRRRPWC (Seq. ID No. 4), CMVHDGADGK (Seq. ID No. 5) and MVHDCADGK (Seq. ID No. 6).

**[0043]** Peptides derived from the carboxyl terminus of urokinase including residues 159-411 were used to generate antibodies specific for the B chain of urokinase also known as the low molecular weight (LMW) form or 33,000 dalton form of urokinase. Six peptides that found use in generating B chain specific antibodies include: YRRHRGGSVTYVC (Seq. ID No. 7), CFIDYPKKEDY (Seq. ID No. 8), SRLNSNTQGEMK (Seq. ID No. 9), SMYNDPQFGTSC (Seq. ID No. 10), LISHRECQPHYYGSEVTTKMLC (Seq. ID No. 11), and SHTKEENGLAL (Seq. ID No. 12).

**[0044]** Other especially useful peptides encompass the peptide bond between lysine 158 and isoleucine 159. These peptides were used to generate antibodies that bind to the urokinase zymogen, also known as the single chain urokinase plasminogen activator (ScuPA). Two peptides that found use in generating urokinase zymogen specific antibodies include: PRFKIIG (Seq. ID No. 13) and RPRFKIIGGE (Seq. ID No. 14). Additional compositions may be considered as immunogens in the preparation of antibodies useful for this invention, for example carbohydrates linked to urokinase or other biochemical or chemical modifications to the urokinase polypeptide chain or peptide sequences homologous to urokinase

sequences. Additionally, combinations of the peptides discussed above, as mixtures, and covalently attached in groups of two or more may be useful.

#### Immunization

[0045] Peptides may be used for immunization directly, after formulation or after conjugation to a carrier molecule. Examples of formulations for peptide immunogen injection include Freund's adjuvant. Carrier molecules that can be used to conjugate to peptide immunogens include bovine serum albumin (BSA), lysozyme, and polysaccharides such as dexetris.

#### Antisera Against Peptides

[0046] Antibodies against the synthetic peptides were raised in mice by two methods (a) immunization by peptide-protein conjugates: Mice were injected and boosted with the peptide-SuBSA conjugate (50 µg per mouse) as an emulsion in complete Freund's adjuvant in the footpads and subcutaneously in the neck. Serial bleedings from 10 days prior, up to 230 days after, the first injection will be collected and studied separately. (b) Immunization with free peptide: In view of the finding (Young et al., *Immunol. Commun.*, 11:9-16 (1982)) that small synthetic peptides (6 residues or larger), when immunized in their free form in complete Freund's adjuvant, will stimulate an *in vivo* antibody response, peptides were injected into mice in their free form (i.e., without coupling to carrier). Each animal received peptide (25-50 µg per mouse) as an emulsion in complete Freund's adjuvant distributed into three sites as above. The animals were boosted with similar doses 3 weeks after the first injection and thereafter monthly. Serial bleedings from 10 days before, up to 200 days after, the initial injection were studied separately.

Monoclonal Antibody Preparation

- [0047] Methods for preparing monoclonal antibodies are well established. See for example, *Monoclonal Antibodies*, eds. Roger H. Kennett, Thomas J. McKearn, Kathleen B. Bechtol, Plenum Press, New York, 1980; *Nature* (1975) 256:495-497.
- [0048] The discovery (Young et al., (1982), *id.*) that synthetic peptides will evoke antibody formation when used as immunogens in their free form (i.e., without coupling to a carrier) was exploited to prepare monoclonal antibodies with preselected submolecular binding specificities to desired protein regions (Schmitz, et al., *Mol. Immunol.*, 19:1699-1702 (1982)). Peptides representing antigenic sites as well as synthetic peptides representing surface regions that are not antigenic when the whole molecule is used as an immunogen have been shown to produce antisera and subsequently monoclonal antibodies of preselected specificities (Schmitz, et al., *Mol. Immunol.*, 20:719-726 (1983); Schmitz, et al., *Mol. Immunol.*, 12:161-175 (1983)).
- [0049] Mice were immunized as above with a given synthetic peptide (25-50  $\mu$ g) in complete Freund's adjuvant and boosted and test-bled at 3 week intervals until high antibody titer was obtained in the test sera. Somatic cell fusions, hybridoma selection, limiting dilution cloning and subcloning, and hybrid cell expansion were performed as described by Schmitz et al. (Schmitz, et al., (1982) *id.*). Expanded subclones were also injected into BALB/cByJ mice ( $2 \times 10^6$  cells/1.0 ml. fresh tissue culture media) that had been primed with pristane (Sigma Chemical Co., St. Louis, MO). Ascites fluid was collected, clarified and stored frozen at -20°C until screened for the presence of hybridoma antibodies.
- [0050] Essentially the peptide immunogen is injected into an animal, then cells

producing a useful composition are immortalized and identified. Additional injections at timed intervals may be used to increase the titer and avidity of the antibodies in the serum. Polyclonal antibodies are prepared by removing serum from these animals. Usually mice are the animal of choice and spleen cells that produce antibodies are immortalized by fusion to immortal mouse myeloma cells, for example SP/2 or 653 cells. Other animals may also be used to generate monoclonal antibodies including chickens or goats or any animal with cells that produce an antibody molecule such that the productive cell can be immortalized. Immortalization of the producing cell is needed to produce a stable starting material for the repeated production of the useful monoclonal antibodies when required. Techniques for cell immortalization include: fusion to immortal cells where the fused cell product continues to express the useful immunoglobulin, or viral infection may be used to immortalize productive cells.

**[0051]** Screening for antibodies with desired specificity can be performed using any of a variety of methods. In a Western blot method, protein of interest is immobilized on a solid surface, such as nitrocellulose, can be incubated with hybridoma supernatants or serum and <sup>125</sup>I-labelled anti-mouse IgG, following removal of the incubation medium and washing of the surface, radioactivity that is bound to the surface is measured. Another method involves addition of the original peptide immunogen to the well of a growing hybridoma culture and looking for blebs or disruptions in the growing cell monolayer. Disruptions are an indication that the hybridoma is producing antibodies against the added immunogen.

**[0052]** Monoclonal antibodies have been prepared which are of use in the present invention. Murine monoclonal antibodies produced from clones A3-E5-A7, A3-E5-B8, A3-E5-G10, A3-E5-H12, A3-E5-B3, A3-E5-D1 were produced against the peptide immunogen PRFKIIG (Seq. ID No. 13). Murine monoclonal antibodies

produced from clones C2-A4-H4, C2-A4-F10, C2-A4-F11 were produced against the peptide immunogen RPRFKIIGGE (Seq. ID No. 14). Murine monoclonal antibodies produced from clones D4-A5-E7, D4-C1-A5, D4-C1-G4, D4-C6-A4, D4-C6-A8, D4-C6-E12, D4, C6-F7 and D4-C6-G12 were directed against the peptide RPRFKIIGGE (Seq. ID No. 14) after it was chemically conjugated to bovine serum albumin.

#### Polyclonal Antibody Preparation

[0053] Essentially any animal that generates an immune response could be used to generate useful polyclonal antibodies for this invention. Animals commonly used include goats, sheep, rats, rabbits, chickens and cows. Particular polyclonal antibody preparations that find use in this invention were prepared from mice and rabbits. Peptides that were of use in generating antibodies against the low molecular weight form of urokinase include YRRHRGGSVTYVC (Seq. ID No. 7), CFIDYPPKKEDY (Seq. ID No. 8), SRLNSNTQGEMK (Seq. ID No. 9), SMYNDPQFGTSC (Seq. ID No. 10), LISHRECQQPHYYGSEVTTKMLC (Seq. ID No. 11), SHTKEENGLAL (Seq. ID No. 12). These antibodies were chemically conjugated to BSA prior to injection. Chemical conjugation is not always required. For example, useful antibodies were also generated against the free peptides CFIDYPPKKEDY (Seq. ID No. 8), SRLNSNTQGEMK (Seq. ID No. 9), and SMYNDPQFGTSC (Seq. ID No. 10) in mice. Other peptides that produced polyclonal antibodies of use in this invention include SNELHQVPSNCD (Seq. ID No. 1), RGKASTDTMGRPCLP (Seq. ID No. 2), CRNPDNRRRP (Seq. ID No. 3), RNPDNRRRPWC (Seq. ID No. 4), CMVHDGADGK (Seq. ID No. 5) and MVHDCADGK (Seq. ID No. 6). Polyclonal antibodies are conveniently prepared by injecting animals with immunogen as described above.



- [0054] Serum is removed from the animal and prepared by conventional techniques as described above. A blood clot is formed spontaneously. The clot is removed and the remaining serum can be purified further. The IgG molecules can be purified away from the serum by affinity chromatography to immobilized protein A or protein G. Proteins that lack an IgG heavy chain constant region will not bind to this affinity resin and will be separated from the desired polyclonal IgG antibodies. Immobilized urokinase may also be used to purify the desired antibody away from other serum proteins. Elution of the bound antibodies can be achieved by using mild denaturants such as 2 molar guanidine hydrochloride or extremes of pH such as a glycine-HCl buffered solution at pH 2.5.
- [0055] The antibody preparations are then screened. Immune IgG preparations or protein A were radiolabeled with  $^{125}\text{I}$  (Amersham Corp., Arlington Heights, IL) using the chloramine-T method (Hunter et al., *Nature*, 194:495-496 (1962)). Unbound  $^{125}\text{I}$  was separated from the radiolabeled sample by gel filtration on Sephadex G-25 (Pharmacia Fine Chemicals, Piscataway, NJ). Protein-associated  $^{125}\text{I}$  was assayed by precipitation with 10% (v/v) trichloroacetic acid.
- [0056] Polyvinylchloride protein assay plates (Costar, Cambridge, MD) were incubated for 3 hours at 37°C with excess (1.5  $\mu\text{g}$  in 50  $\mu\text{g}$  of PBS/well) test and control antigens, washed extensively with PBS, and blocked with 1% BSA in PBS (100  $\mu\text{l}$ /well) for 1 hr at 37°C to prevent non-specific binding of subsequent reagents. After washing, the plates were used for binding antibody.
- [0057] Sera, culture supernatants and clarified ascites fluids were screened for anti-UK antibodies using a solid phase RIA described by Sakata and Atassi (*Mol. Immunol.*, 18:961-967 (1981) as modified by Schmitz et al. ((1982) *id.*). This assay was also used to determine antibody binding specificities to peptides. Briefly, RIA plates that had been coated with the appropriate test antigens (various

forms of UK or peptide conjugates) were incubated for 3 hours at 37°C with an antibody preparation (50 µl/well) appropriately prediluted in PBS-BSA so as to maximize specific binding. The plates were subsequently washed with PBS and amplified with excess (1:1000 dilution of the stock reagent in PBS-BSA) rabbit anti-mouse IgG + IgM antisera (Litton Bionetics, Kensington, MS) for 2 hrs at 37°C. After washing, the plates were developed with excess ( $2.0 \times 10^5$  cpm in 50 µl PBS-BSA/well)  $^{125}$ I-labelled protein A for 2 hrs at room temperature, washed, and then separated into individual wells that are counted in a gamma counter (Beckman Instruments, Inc., Irvine, CA). Results are corrected for nonspecific (0.1-2%) binding detected in control wells not coated with test antigen but blocked with BSA.

[0058] When peptide-SuBSA conjugates are the immunizing antigens then antibody responses were analyzed by peptide-lysozyme (not succinyl lysozyme) conjugates since lysozyme and SuBSA do not cross-react immunochemically. Correction for non-specific binding is derived from binding to lysozyme control. When the free peptides were the antigens then antibodies were analyzed on peptide-SuBSA conjugates and employing lysozyme and SuBSA as controls.

[0059] One of the novel aspects of these antibodies is their specificity for binding only to urokinase and not to other related proteins such as trypsin or TPA.

#### Other Methods of Using the Compositions

##### Analysis of Separated Kidney Cells

[0060] Cell separation in microgravity is one of the significant new applications of space technology to basic and clinical medical research. Separations are more effective in microgravity conditions due to the lack of sedimentation and absence of other limitations imposed by physical forces. Kidney cells have been used as a

model for their separation by continuous flow electrophoresis and other methods under microgravity conditions. After the cells are separated, it is important to determine their purity, viability and physiological performance. The compositions and methods provided in this invention have been used to determine these parameters.

#### Cancer Diagnosis and Monitoring

[0061] This invention provides methods and compositions that are useful in the diagnosis and monitoring of tumors, which produce urokinase (such as breast, prostate and adenosarcoma). Analysis can be performed on biological fluids such as blood serum, urine, saliva and semen as a way of detecting urokinase and identifying cancer development or progress. The cells are removed from the serum by centrifugation and the soluble protein fraction assayed by a sensitive ELISA or RIA technique as previously described.

[0062] Immunohistochemical methods can be used to identify the presence of urokinase in tumor biopsy specimens. For example, the tumor tissue can be quick frozen, sections can be cut and fixed with acetone on glass slides, the sections can be treated with normal serum to reduce non-specific anti-urokinase antibody binding. The antigen can be detected by exposure to unlabeled primary antibody and a labeled second antibody. Cells from fine needle biopsies may be similarly tested for uPA.

Purification of Individual Forms of Urokinase

[0063] Antibodies are relatively stable macromolecules they can be immobilized on solid supports to produce useful chromatography media. The unique specificity of an antibody makes the chromatography media extremely specific for particular molecules and facilitates the purification of these molecules from complex mixtures. A variety of methodologies can be used to immobilize antibodies on solid supports. Antibodies useful to this invention can be used in this way to provide chromatography resin that will allow the purification of each form of urokinase.

[0064] The following examples are offered as illustrations and are not meant to limit the scope of the invention.

**EXAMPLES****Example I: Synthesis and Purification of Peptide Immunogens**

[0065] A total of 14 peptides corresponding to different areas of the urokinase molecule have been synthesized, purified and characterized. A general discussion of the chemical synthesis of peptides, their purification and characterization can be found in the Description of Preferred Embodiments above.

**Peptides Synthesized for Antibody Development**

- |      |         |  |
|------|---------|--|
| (1)  | 1-12    | SNELHQVPSNCD (Seq. ID No. 1)             |
| (2)  | 59-73   | RGKASTDTMGRPCLP (Seq. ID No. 2)          |
| (3)  | 102-111 | CRNPDNRRRP (Seq. ID No. 3)               |
| (4)  | 103-113 | RNPDNRRRPWC (Seq. ID No. 4)              |
| (5)  | 126-135 | CMVHDGADGK (Seq. ID No. 5)               |
| (6)  | 127-135 | MVHDCADGK (Seq. ID No. 6)                |
| (7)  | 177-189 | YRRHRGGSVTYVC (Seq. ID No. 7)            |
| (8)  | 205-215 | CFIDYPPKKEDY (Seq. ID No. 8)             |
| (9)  | 222-233 | SRLNSNTQGEMK (Seq. ID No. 9)             |
| (10) | 319-341 | SMYNDPQFGTSC (Seq. ID No. 10)            |
| (11) | 282-293 | LISHRECQQPHYYGSEVTTKMLC (Seq. ID No. 11) |
| (12) | 401-411 | SHTKEENGLAL (Seq. ID No. 12)             |
| (13) | 155-161 | PRFKIIG (Seq. ID No. 13)                 |
| (14) | 154-163 | RPRFKIIGGE (Seq. ID No. 14)              |

**Example II: Preparation of Polyclonal Antibodies**

[0066] A general description of the method used to prepare polyclonal antibodies can be found in Description of the Preferred Embodiments above. Antisera obtained in this way binds very strongly to the immunizing peptides as shown by solid phase RIA.

[0067] The polyclonal antisera were labelled with  $^{125}\text{I}$  as described above and their binding to the peptide-lysozyme conjugate and to the low molecular weight urokinase was determined. Results represent pools of antisera from 3 mice each and varied  $\pm 1.7\%$  or less.

**Table 1. Antibodies against UK Peptides of the B chain bind to both correlative peptide-lysozyme conjugates and to low molecular weight urokinase, with similar affinity.**

Antigen	$^{125}\text{I}$ -Antibodies bound ( $\Delta$ cpm) to	
	Peptide-lysozyme conjugate	Low M.W. U.K.
peptide 177-189 (Seq. ID No. 7)	43,940	45,570
peptide 205-215 (Seq. ID No. 8)	42,490	44,310
peptide 222-233 (Seq. ID No. 9)	45,570	47,520
peptide 319-330 (Seq. ID No. 10)	49,540	51,090
peptide 319-341 (Seq. ID No. 11)	42,590	43,510
peptide 401-411 (Seq. ID No. 12)	41,060	42,520

[0068] It was similarly shown that certain of the C-terminal domain peptides demonstrated binding to all forms (LMW and HMW) of urokinase with affinities comparable to the binding exhibited by the correlative free peptide.

**Table 2. Antisera (polyclonal) Against Peptides in the common C-terminal Domain Bind All Forms of Urokinase.****SPECIFIC BINDING (CPM)**

<u>Antigen</u>	<u>Mouse No.</u>	<u>Binding to Free Peptide</u>	<u>Binding to LMW</u>	<u>Urokinase HMW</u>
Peptide 205-215 (Seq. ID No. 8)	1	50033	30761	44755
	2	24190	10563	17550
	3	40467	23841	36670
Peptide 222-233 (Seq. ID No. 9)	1	36101	17354	26819
	2	52730	19853	27951
	3	25601	10635	17033
Peptide 319-330 (Seq. ID No. 10)	1	14863	10923	11699
	2	24374	9617	16822
	3	9183	5168	6833

**[0069]** The same test were run on N-terminal domain peptides with similar results against whole urokinase.

Table 3. Antisera (polyclonal) Binding to Urokinase Peptides and to Urokinase.

**SPECIFIC BINDING**

<u>Antigen</u>	<u>Mouse No.</u>	<u>Binding to Free Peptide</u>	<u>Binding to Urokinase</u>
Peptide 1-12 (Seq. ID No. 1)	1	30740	8726
	2	36318	3644
	3	27112	3554
Peptide 59-73 (Seq. ID No. 2)	1	33478	3554
	2	19784	3294
	3	25172	3358
Peptide 102-111 (Seq. ID No. 3)	1	25992	10414
	2	21596	6482
	3	36336	8988
Peptide 103-113 (Seq. ID No. 4)	1	18876	7220
	2	10964	8708
	3	12656	9788
Peptide 126-135 (Seq. ID No. 5)	1	35798	3784
	2	15668	2016
	3	21434	3460

**Example III: Preparation of Monoclonal Antibodies**

[0070] A general description of the method used for the preparation of monoclonal antibodies can be found in Description of Preferred Embodiments. Monoclonal antibodies were prepared with the capacity to discriminate between the three common molecular forms of urokinase. Three groups of antibodies were



developed by the use of synthetic peptides which represent the areas described above. The first group of peptides consisted of areas present both in prourokinase and high molecular weight urokinase, but absent in the low molecular weight form. These six peptides were synthesized, cleaved, purified and injected into BALB/c mice for the preparation of peptide-specific antisera. Spleen cells from these mice were fused with SP/2 myeloma cells in order to prepare monoclonal antibodies. In addition, six peptides representing areas present in all three molecular forms were made. Monoclonal antibodies raised against these peptides bound all forms of urokinase but not other serine proteases.

[0071] For the differential recognition of prourokinase from high molecular weight urokinase, specific monoclonal antibodies, binding to the area of cleavage, were developed. Two peptides were constructed which represent the area around the single peptide bond that is cleaved in the activation process (i.e., conversion of the zymogen into the active enzyme). Monoclonal antibodies specific for this area only bound to urokinase zymogen but not high molecular weight urokinase.

[0072] After injection of the peptides into outbred mice, and fusion of their spleen cells, 17 monoclonal antibodies were selected for their specificity to one of the peptides (peptide 154-163, Seq. ID No. 14). A table showing the results of this screening is shown at Table 4, below.

[0073] These monoclonals were tested for binding to the different forms of urokinase by the western blot method. The nitrocellulose paper was incubated with cell culture supernatants and <sup>125</sup>I-labelled anti-mouse IgG. Three of the monoclonals showed specific binding to single chain urokinase. The clones involved were expanded and the monoclonal antibodies purified for further characterization and testing.

[0074] Several monoclonal antibodies were also made in mice and were then tested for specificity to ScuPA, the A-chain, B-chain the cleavage bonds at Lys-158—Ile-

159, the loss of Lys-158 and the cleavage bond at Lys-135—Lys136 during activation. Additional murine monoclonal antibodies and rabbit anti-urokinase monoclonal antibodies were made and tested.

**Table 4. Antibody Titer of Urokinase Monoclonals (cpm/min of RIA)**

<u>Clone No.</u>	<u>BSA-Urokinase Peptide</u>	<u>BSA-Non Sense Peptide</u>	<u>BSA Only Blank</u>
A3-E5-A7	71,813	803	901
A3-E5-B8	47,816	1,321	876
A3-E5-G10	58,503	1,454	1,311
A3-E5-H12	85,240	1,138	1,158
A3-E5-B3	59,674	871	821
A3-E5-D1	36,729	934	877
C2-A4-H4	7,633	1,855	928
C2-A4-F10	40,659	3,053	1,448
C2-A4-F11	91,196	6,987	1,726
D4-A5-E7	4,962	1,393	1,430
D4-C1-A5	25,631	1,579	1,355
D4-C1-G4	90,645	3,410	1,497
D4-C6-A4	72,303	1,101	1,045
D4-C6-A8	43,309	946	2,171
D4-C6-E12	82,320	1,405	1,055
D4-C6-F7	68,513	1,280	868
D4-C6-G12	43,013	1,109	781

Note: Group A; mice were immunized with free peptide 153-162 (Seq. ID No. 15)

Group C; mice were immunized with free peptide 154-163 (Seq. ID No. 14)

Group D; mice were immunized with BSA-conjugated peptide 154-163  
(Seq. ID No. 14)

**Example IV: Comparison of Prior Art Test for Urokinase with That of Invention**

[0075] Comparison of three assay methods for measuring urokinase in supernatants from kidney cell cultures were made in order to test the efficacy of the methods

and compositions of the invention with those of the prior art. Cell cultures of HEK cells were prepared as described previously above and supernatants of the cultures were isolated. The supernatants were divided into aliquots for standardization runs, and runs using one of the three test methods.

[0076] The S-244 method used a chromatogenic substrate to measure active urokinase. The method used was that of Claeson, G. et al., "Methods for Determination of Prekallikrein in plasma, Glandular Kallikrein, and Urokinase," *Haemostasis* 7:76-87 (1978), incorporated specifically herein to the extent that it provides materials and methods not otherwise provided herein. The microclot lysis assay ("MCLA") assay measures the optical density of the clots being lysed by the urokinase to determine the half-lysis time for standard concentrations. The method used was that of Lewis et al., "A Miniaturized Fibrinolytic Assay for Plasminogen Activators," *Thrombosis Res.* 64:223-234 (1991), incorporated specifically herein to the extent that it provides materials and methods not otherwise provided herein.

[0077] The molecular-specific antibody ("MSA") assay of the invention has been previously described. It is typically run as an RIA ("MSA(RIA)") method (although it can also be an ELISA assay or other non-radioactive assay) using antibodies which are specific for the intact region of the active site (158-159) on the urokinase molecule. The RIA methods utilized were similar to those previously described above for the screening of supernatants for anti-UK antibodies using solid-phase assays (Sakata and Atassi 1981, as modified by Schmitz et al. 1982). While it is possible to label either the supernatant proteins to be tested or to label the antibodies, the inventors typically label the antibodies as previously described.

[0078] All assays were calibrated against an international reference standard for urokinase (IRP-UK; National Biologics Standards Board, London). A typical

standard curve for the RIA assay (specifically, those used to standardize values observed for the Run 4 Day 12 supernatants below) is shown in Fig. 1. That standard curve plots the results from three separate standard assays as net counts per minute versus microgram concentrations of urokinase in 50 microliter volumes.

[0079] Fig. 2 is a graph comparing the amount of scuPA in the Run 4 Day 12 supernatants using the differential S-2444 assay, the MCLA assay, and the MSA(RIA) method of the invention. The amount of scuPA contained in the supernatant culture medium harvested from four (4) day old cultures of human kidney cells, which cells had been separated by continuous flow electrophoresis into twenty-four (24) subpopulations, were determined. The estimates of scuPA found by the S-2444 assay were obtained by using two different aliquots of the supernatant and subtracting the levels of uPA activity found in the S-2444 assay (direct) and the levels found in the S-2444 assay following pre-incubation with 15 nM plasmin (2 hours at 37°C), (which converts all of the scuPA to active two-chain uPA). This shows that the pre-incubation with plasma misses a significant amount of scuPA which is actually present in the culture medium. Conversely, the MSA is sensitive enough to measure it directly. Similarly, the MCLA assays routinely detected far less of the scuPA than did the MSA(RIA) methods of the invention.

[0080] Typically, the MCLA assay is the most comparable to the MSA(RIA) assay. Thus, for instance as can be seen in supernatant fraction 24, the two values are almost identical. Conversely, both the MCLA and the MSA(RIA) methods routinely detect substantially higher levels of scuPA than does the S-2444 assay. This was not always the case, however, since, for instance, in supernatant fractions 2-4, the S-2444 assay detected slightly more scuPA than did the MCLA

assay. However, in all instances, the MSA(RIA) assays of the invention were able to detect substantially higher levels of scuPA than the two prior art approaches.

\* \* \* \* \*

[0081] The present invention has been described in terms of particular embodiments found or proposed to comprise preferred modes for the practice of the invention. It will be appreciated by those of skill in the art that, in light of the present disclosure, numerous modifications and changes can be made in the particular embodiments exemplified without departing from the intended scope of the invention. For example, the antibody preparations of the invention would be useful for whole body imaging techniques to look for tumors with elevated local concentrations of urokinase. Immunotoxins based on these antibody preparations would likewise be useful. Fragments such as Fab or Fv derived from the invention antibodies will find use in certain instances, as well. The peptides might also be used as cancer vaccines to activate the human immune system against tumors expressing urokinase. Additionally, immunoaffinity columns may be constructed with the discriminating antibodies of the invention to quantitate the different forms of urokinase in a sample (e.g., *See, Sibley, et al., Biotechniques* Vol. 10, 1993; Universal Sensors, Metairie, LA). All such modifications are intended to be included within the scope of the appended claims.

SEQUENCE LISTING

(1) GENERAL INFORMATION

- (i) APPLICANT: Atassi, M. Z.  
Morrison, D.R.
- (ii) TITLE OF INVENTION: Molecular-Specific Urokinase  
Antibodies
- (iii) NUMBER OF SEQUENCES: 16
- (iv) CORRESPONDENCE ADDRESS:
  - (A) ADDRESSEE: NASA Johnson Space Center
  - (B) Mail Code: HA/Office of Patent Counsel
  - (C) STREET: 2101 NASA Road 1
  - (D) CITY: Houston
  - (E) STATE: Texas
  - (F) COUNTRY: USA
  - (G) ZIP: 77546
- (v) COMPUTER READABLE FORM:
  - (A) MEDIUM TYPE: Floppy disk
  - (B) COMPUTER: IBM PC compatible
  - (C) OPERATING SYSTEM: PC-DOS/MS-DOS
  - (D) SOFTWARE: PatentIn Release #1.0,  
Version #1.25
- (vi) CURRENT APPLICATION DATA:
  - (A) APPLICATION NUMBER:
  - (B) FILING DATE:
  - (C) CLASSIFICATION:
- (viii) ATTORNEY/AGENT INFORMATION:
  - (A) NAME: Cate, James M.
- (ix) TELECOMMUNICATION INFORMATION:
  - (A) TELEPHONE: 281-483-1001
  - (B) TELEFAX: 281-244-8452

## (2) INFORMATION FOR SEQ ID NO:1:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 12 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

## (ii) MOLECULE TYPE: peptide

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO:1:

Ser Asn Glu Leu His Gln Val Pro Ser Asn Cys Asp  
1 5 10

## (2) INFORMATION FOR SEQ ID NO:2:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 15 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

## (ii) MOLECULE TYPE: peptide

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO:2:

Arg Gly Lys Ala Ser Thr Asp Thr Met Gly Arg Pro Cys Leu Pro  
1 5 10 15

## (2) INFORMATION FOR SEQ ID NO:3:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 10 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

## (ii) MOLECULE TYPE: peptide

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO:3:

Cys Arg Asn Pro Asp Asn Arg Arg Arg Pro  
1 5 10

## (2) INFORMATION FOR SEQ ID NO:4:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 11 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

## (ii) MOLECULE TYPE: peptide

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO:4:

Arg Asn Pro Asp Asn Arg Arg Arg Pro Trp Cys  
1 5 10

## (2) INFORMATION FOR SEQ ID NO:5:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 10 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: circular

## (ii) MOLECULE TYPE: peptide

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO:5:

Cys Met Val His Asp Gly Ala Asp Gly Lys  
1 5 10

## (2) INFORMATION FOR SEQ ID NO:6:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 9 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: circular

## (ii) MOLECULE TYPE: peptide

## (xi) SEQUENCE DESCRIPTION: SEQ ID NO:6:

Met Val His Asp Cys Ala Asp Gly Lys  
1 5



(2) INFORMATION FOR SEQ ID NO:7:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 13 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:7:

Tyr	Arg	Arg	His	Arg	Gly	Gly	Ser	Val	Thr	Tyr	Val	Cys
1				5					10			

(2) INFORMATION FOR SEQ ID NO:8:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 11 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:8:

Cys	Phe	Ile	Asp	Tyr	Pro	Lys	Lys	Glu	Asp	Tyr
1				5					10	

(2) INFORMATION FOR SEQ ID NO:9:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 12 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:9:

Ser	Arg	Leu	Asn	Ser	Asn	Thr	Gln	Gly	Glu	Met	Lys
1				5					10		

(2) INFORMATION FOR SEQ ID NO:10:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 12 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:10:

Ser Met Tyr Asn Asp Pro Gln Phe Gly Thr Ser Cys  
1 5 10

(2) INFORMATION FOR SEQ ID NO:11:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 23 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:11:

Leu Ile Ser His Arg Glu Cys Gln Gln Pro His Tyr Tyr Gly Ser Glu  
1 5 10 15  
Val Thr Thr Lys Met Leu Cys  
20

(2) INFORMATION FOR SEQ ID NO:12:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 11 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:12:

Ser	His	Thr	Lys	Glu	Glu	Asn	Gly	Leu	Ala	Leu
1				5					10	

(2) INFORMATION FOR SEQ ID NO:13:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 7 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:13:

Pro	Arg	Phe	Lys	Ile	Ile	Gly
1				5		

(2) INFORMATION FOR SEQ ID NO:14:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 10 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:14:

Arg	Pro	Arg	Phe	Lys	Ile	Ile	Gly	Gly	Glu
1				5					10

(2) INFORMATION FOR SEQ ID NO:15

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 10 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:15

Leu	Arg	Pro	Arg	Phe	Lys	Ile	Ile	Gly	Gly
1				5				10	

(2) INFORMATION FOR SEQ ID NO:16:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 411 amino acids

(B) TYPE: amino acid

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:16:

Ser	Asn	Glu	Leu	His	Gln	Val	Pro	Ser	Asn	Cys	Asp	Cys	Leu	Asn	Gly
1				5				10						15	
Gly	Thr	Cys	Val	Ser	Asn	Lys	Tyr	Phe	Ser	Asn	Ile	His	Trp	Cys	Asn
			20					25					30		
Cys	Pro	Lys	Lys	Phe	Gly	Gly	Gln	His	Cys	Glu	Ile	Asp	Lys	Ser	Lys
		35					40					45			
Thr	Cys	Tyr	Glu	Gly	Asn	Gly	His	Phe	Tyr	Arg	Gly	Lys	Ala	Ser	Thr
	50					55					60				
Asp	Thr	Met	Gly	Arg	Pro	Cys	Leu	Pro	Trp	Asn	Ser	Ala	Thr	Val	Leu
65					70					75					80
Gln	Gln	Thr	Tyr	His	Ala	His	Arg	Ser	Asp	Ala	Leu	Gln	Leu	Gly	Leu
				85					90					95	
Gly	Lys	His	Asn	Tyr	Cys	Arg	Asn	Pro	Asp	Asn	Arg	Arg	Arg	Pro	Trp
			100					105						110	
Cys	Tyr	Val	Gln	Val	Gly	Leu	Lys	Pro	Leu	Val	Gln	Glu	Cys	Met	Val
		115					120					125			
His	Asp	Cys	Ala	Asp	Gly	Lys	Lys	Pro	Ser	Ser	Pro	Pro	Glu	Glu	Leu
	130					135					140				

**MSC-21947-1-CU**

**Express mailing certificate No. EU657647915US**  
**Patent Application**

Lys 145	Phe	Gln	Cys	Gly	Gln 150	Lys	Thr	Leu	Arg	Pro 155	Arg	Phe	Lys	Ile	Ile 160
Gly	Gly	Glu	Phe	Thr 165	Thr	Ile	Glu	Asn	Gln 170	Pro	Trp	Phe	Ala	Ala	Ile 175
Tyr	Arg	Arg	His 180	Arg	Gly	Gly	Ser	Val 185	Thr	Tyr	Val	Cys	Gly	Gly	Ser 190
Leu	Ile	Ser 195	Pro	Cys	Trp	Val	Ile	Ser 200	Ala	Thr	His	Cys	Phe	Ile	Asp 205
Tyr	Pro 210	Lys	Lys	Glu	Asp	Tyr	Ile	Val	Tyr	Leu	Gly 220	Arg	Ser	Arg	Leu
Asn 225	Ser	Asn	Thr	Gln	Gly 230	Glu	Met	Lys	Phe	Glu 235	Val	Glu	Asn	Leu	Ile 240
Leu	His	Lys	Asp	Tyr 245	Ser	Ala	Asp	Thr	Leu	Ala 250	His	His	Asn	Asp	Ile 255
Ala	Leu	Leu	Lys 260	Ile	Arg	Ser	Lys	Glu 265	Gly	Arg	Cys	Ala	Gln	Pro	Ser 270
Arg	Thr 275	Ile	Gln	Thr	Ile	Cys	Leu	Pro 280	Ser	Met	Tyr	Asn	Asp	Pro	Gln 285
Phe 290	Gly	Thr	Ser	Cys	Glu	Ile 295	Thr	Gly	Phe	Gly	Lys 300	Glu	Asn	Ser	Thr
Asp 305	Tyr	Leu	Tyr	Pro	Glu 310	Gln	Leu	Lys	Met	Thr 315	Val	Val	Lys	Leu	Ile 320
Ser	His	Arg	Glu	Cys 325	Gln	Gln	Pro	His	Tyr 330	Tyr	Gly	Ser	Glu	Val	Thr 335
Thr	Lys	Met	Leu 340	Cys	Ala	Ala	Asp	Pro 345	Gln	Trp	Lys	Thr	Asp	Ser	Cys 350
Gln	Gly	Asp 355	Ser	Gly	Gly	Pro	Leu	Val 360	Cys	Ser	Leu	Gln	Gly	Arg	Met 365
Thr 370	Leu	Thr	Gly	Ile	Val	Ser	Trp	Gly 375	Arg	Gly	Cys	Ala	Leu	Lys	Asp 380

MSC-21947-1-CU

Express mailing certificate No. EU657647915US  
Patent Application

Lys Pro Gly Val Tyr Thr Arg Val Ser His Phe Leu Pro Trp Ile Arg  
385 390 395 400

Ser His Thr Lys Glu Glu Asn Gly Leu Ala Leu  
405 410

(2) INFORMATION FOR SEQ ID NO:17

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 8 amino acids
- (B) TYPE: amino acid
- (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: peptide

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:17

Ala Asp Gly Lys Lys Pro Ser Ser  
1 5